

Table 1. Suppl. List of different adapters and primers used in methylation-sensitive amplification polymorphism analysis.

Adapter/primer	Nucleotide sequences
<i>EcoRI</i> adapter	5'-CTC GTAGACTGCGTACC-3' 3'-CATCTGACGCATGGTTAA-5'
<i>MspI-HpaII</i> adapters	5'-GA TCATGAGTCCTGCT-3' 3'-AGTACTCAGGACGAGC-5'
Pre-amplification primers	
<i>EcoRI</i>	5'-GACTGCGTACCAATTCA-3'
<i>MspI-HpaII</i>	5'-ATCATGAGTCCTGCTCGG-3'
Selective amplification primers	
E + 2 TT	5'-GACTGCGTACCAATTCATT-3'
TG	5'-GACTGCGTACCAATTCATG-3'
MH + 3 AAC	5'-ATCATGAGTCCTGCTCGGAAC-3'
ATG	5'-ATCATGAGTCCTGCTCGGATG-3'
AAG	5'-ATCATGAGTCCTGCTCGGAAG-3'
ACA	5'-ATCATGAGTCCTGCTCGGACA-3'